

EMULADOR DE ESTERILIZACIÓN PARA VAPOR TIPO 6 - Sv = 134 °C - 3.5 min.

REF. 4500014009

Tiras indicadoras de verificación de ciclo tipo 6 para procesos de esterilización para vapor.

Resultado inmediato

Precisión de acuerdo con la norma ISO 11140 para indicadores tipo 6.

Valores establecidos mínimos (Sv): 134 °C - 3.5 minutos


Libre de plomo y de metales tóxicos

Color inicial: morado

Color final: verde

Presentación: Tiras de papel bilaminado de 75 x 21.5 mm con patrón de viraje en cada tira. Envase x 200 unidades.

Marca Kims - Industria Argentina



The product packaging for the Type 6 sterilization indicator strip is shown. It features a large blue number "6" on the left side. The text on the packaging includes: "Emulador de esterilización TIPO 6", "TYPE 6 Sterilization Emulator", "Emulador para esterilização TIPO 6", "VAPOR [STEAM]", "Contenido - Contents - Contém 200 unidades / units", "Made in Argentina", and "Valores establecidos Standard values 134 °C - 3.5 min". The KIMS logo is also present. A circular logo with "Pb" and "Libre de plomo" (Lead free) is located on the left side. In the center, there is a diagram showing a purple circle changing to a green circle, with an arrow pointing from the purple circle to the green circle. Below the diagram, the text reads: "TURNS FROM PURPLE TO GREEN", "VIRA DE MORADO A VERDE", and "MUDANÇA DE COR: PURPURA PARA VERDE".